



CALCULICLAST



Extra power for stone fragmentation precise & minimum forward movement of the probe minimizes stone flying.

Power of the blast adjustable. Four preset frequencies of the blast.

Special accessories & probes for fragmentation of renal & bladder stones. Adaptors to connect the machine to central connection.

New Suction Cannula with Pinch Valve

INPUT VOLTAGE 230 Volts. **OUTPUT PRESSURE** 4kg/cm2 **AIRFLOW**

TECHNICAL SPECIFICATIONS

Control Unit

Supply Voltage - 230 Volts A.c. Working pressure - 1.8 to 4 kg/cm2. Maximum input pressure - Max. 6.0 Kg/sq.cm. Pressure Source - Pressurised Oxygen / Compressed Air Hand Piece

Dimensions - Dia. 24 Mm X 180 Mm. Weight - 0.150 Gms.

Probes

0.8 mm .1 mm .1.2mm for ureternoscopy 1.5 mm Bladder Stones. 2 mm renal stones.



CALCULICLAST WORKSTATION

CALCULICLAST with oil less, noiseless compressor & powder coated trolley.

FEATURES

30 liter/min.

- Oil less, noiseless compressor
- Electrical cut off switch
- Mechanical safety valve
- Powder coated trolley with wheels

Maintenance Free Design.

INPUT VOLTAGE

230 Volts,

OUTPUT PRESSURE

4kg/cm2

AIRFLOW

30 liter/min.



CALCULICLAST

ASTRA 30 - HO:YAG LASER Lithotripter

URODYNAMICS

UROFLOWMETER

WORK STATION

HYDROMAT

WIRELESS UROFLOWMETER

SELFCLEANING UROFLOWMETRY



STATUS MEDICAL EQUIPMENTS INDIA PVT LTD

M-9/A, Additional M.I.D.C.

Satara 415004,

Maharashtra India

Phone 91-02162-240037 e-mail info@statmedi.com

www.statmedi.com





ISO 13485:2016 Certified Company